

Identification of *Alternaria alternata* in combination with analysis of the G143A substitution

June 2014

This method was developed by Byron Vega and Megan Dewdney (published in 2014: Distribution of Qol resistance in populations of Tangerine-infecting *Alternaria alternata* in Florida. Plant Disease 98, 67-76.). The PCR-conditions below and the digestion of PCR-products are adapted to suit the reaction enzymes used.

The advantage with this method is that it is species specific for *A. Alternata* and that the PCR product can be cut with a restriction enzyme to reveal the amino acid at position 143 in the gene encoding cytochrome *b*, associated with loss of sensitivity against strobilurins.

DNA extraction

- Cut out a leaf disc containing one single lesion (*ca.* 5 mm \varnothing), which include both green and necrotic tissue.
- Washed it twice in sterile distilled water and dried with paper towels.
- Homogenised using five glass beads (3mm) in a 2 ml micro centrifuge tube and shake at 5500 repeats per minute during 30 s (*e.g.* Precellys® 24 Bead Mill Homogenizer, Bertin Technologies).
- Add the extraction buffer. (800 mL of 3% CTAB-EDTA, pH 8 or the amount recommended if using a kit). Shake again as above.
- Follow the CTAB-EDTA protocol or the kit protocol onwards.
- Re-suspend the DNA pellet in 30 μ L distilled water if CTAB-protocol or as written in the kit-protocol.
- Dilute the samples to approx. 1 ng DNA μ l⁻¹.

PCR protocol

Forward primer: Cytb2f 5'-CTA TGG ATC TTA CAG AGC AC-3'

Reverse primer: DTRcytb2-INTr 5'-GTA TGT AAC CGT CTC CGT C-3'

PCR solution of 20 μ l contains 10 μ l DNA, 2.75 mM MgCl₂ (final concentration), 0.2 mM dNTP, 0.02 μ M of each primer, 0.05 U μ l⁻¹ of DreamTaq® DNA Polymerase (Fermentas International Ink, Canada) and 10X DreamTaq™ Green Buffer.

The amplification product is 377 bp, profile I according to Vega and Dewdney (2014). They have found a longer product, 1,564 bp, for profile II, but this has not been observed in Sweden so far.

Master mix 10 μ l + 10 μ l DNA:	H ₂ O	4,5 μ l
	Buffer (10x)	2 μ l
	dNTP (2 mM)	2 μ l
	MgCl (25 mM)	0.6 μ l
	Primer F (10 nM)	0.4 μ l
	Primer R (10 nM)	0.4 μ l
	Taq (5 U μ l ⁻¹)	0.1 μ l

PCR condition: 95°C for 2 min, 40 cycles of 30 s at 95°C, 30 s at 59°C and 30 s at 72°C, followed by a 5 min extension at 72°C. Separate the products with electrophoresis on 1% agarose gel stained with suitable dye and visually analyse as under UV-light.

The presence of the G143A substitution is tested using Cleaved Amplified Polymorphic Sequences (CAPS) where the PCR product is digested using the restriction enzyme *SatI* (*Fnu4HI*) (Thermo Scientific Inc.). Incubate the reaction volume of 10 µl containing 5 U of enzyme solution and 5 µl of PCR product solution in a 37°C water bath for two hours. Separate the PCR products in a 1.5% agarose gel.

If the G143A is present in the sample, the length of the digested fragments will be 123 bp and 254 bp for profile I and 254 bp and 1,310 bp for profile II. The fragments of 377 and 1,564 bp remain undigested in wild-type isolates for profile I and profile II respectively.

Sent in by Eva Edin, Dept. Forest Mycology & Plant Pathology, Swedish University of Agricultural Sciences, P.O. Box 7026, SE 75007 Uppsala Sweden. eva.edin@slu.se